


AddOn[®] hydrophilic · EDOF


Technical specifications

Refractive Index	1.46	Lens type	One piece sulcus lens
Abbe No.	58	Material	Biocompatible hydrophilic acrylic (25%)
		Haptic	4 Flex-Haptic
		Optic	spheric
		Optic form	convex-concave
		Filter	UV filter
		Sterilisation	Steam
		Injection system	1stInject 2.0HB, 1stInject 2.4HB



■ Basis function

Optical Power




■ EDOF function

low addition

AddOn[®] hydrophilic · EDOF

Article No.	Powers	Steps	Addition
A4DW0H	0.0 D		low addition
A4EW0H	-5.0 D to -0.5 D	0.25 D	low addition
A4EW0H	+0.5 D to +5.0 D	0.25 D	low addition



For the latest version,
scan the QR code.
1stq.de/addon1069

**Do you have any
questions?**
export@1stq.eu



1stQ GmbH
Konrad-Zuse-Ring 23
68163 Mannheim, Germany
www.1stq.de/en