


AddOn[®] hydrophilic · toric



Technical specifications

Refractive Index 1.46

Abbe No. 58

Lens type One piece sulcus lens

Material Biocompatible hydrophilic acrylic (25%)

Haptic 4 Flex-Haptic

Optic spheric

Optic form convex-concave

Filter UV filter

Sterilisation Steam

Injection system 1stInject 2.0HB to Cyl. 4.50 D
1stInject 2.4HB

■ Basis function

Optical Power (spherical equivalent)



■ Cylinder function

Cylindrical Power



AddOn[®] hydrophilic · toric

Article No.	SEQ	Steps	Cylinder	Steps
A4TW0T	0.0 D		1.00 D	0.50 D
A4TW0T	0.0 D		1.50 D to 4.50 D	0.75 D
A4TW00	0.0 D		5.25 D to 8.25 D	0.75 D
A4TW00	0.0 D		9.00 D to 11.00 D	1.00 D
A4FW0T	-10.0 D to +10.0 D (not 0.0 D)	0.25 D	1.00 D	0.50 D
A4FW0T	-10.0 D to +10.0 D (not 0.0 D)	0.25 D	1.50 D to 4.50 D	0.75 D
A4FW00	-10.0 D to +10.0 D (not 0.0 D)	0.25 D	5.25 D to 8.25 D	0.75 D
A4FW00	-10.0 D to +10.0 D (not 0.0 D)	0.25 D	9.00 D to 10.00 D	1.00 D
A4FW00	-3.0 D to +8.0 D	0.25 D	11.00 D	



For the latest version,
scan the QR code.
1stq.de/addon1052

**Do you have any
questions?**
export@1stq.eu



1stQ GmbH
Konrad-Zuse-Ring 23
68163 Mannheim, Germany
www.1stq.de/en