



Technical specifications

Haigis-constants*	a0 = -0.7633 a1 = 0.1391 a2 = 0.2244
A-constant*	118.895 (SRK/T)
Refractive Index	1.5
Abbe No.	50

*IOL Con 12.2021

Lens type	One piece capsular bag IOL
Material	Biocompatible hydrophobic acrylic (<0.5%)
Haptic	C Haptic, 0° angulation
Optic	Aspheric optic, 360° cell barrier
Filter	UV filter, blue light filter (optional)
Sterilisation	Ethylene oxide
Injection system	1stInject 2.0HB, 1stInject 2.4HB / Medicel Accuject

Basis function

Optical Power



Basis V	Art. No. - clear	Art. No. - yellow	Powers	Steps
non-preloaded	B5AD00	B5ADY0	0.0 D to +9.0 D	1.0 D
non-preloaded	B5AD00	B5ADY0	+10.0 D to +30.0 D	0.5 D
preloaded	B5AB00	B5ABY0	0.0 D to +9.0 D	1.0 D
preloaded	B5AB00	B5ABY0	+10.0 D to +30.0 D	0.5 D



For the latest version,
scan the QR code.
1stq.de/basisv1065

**Do you have any
questions?**
export@1stq.eu



1stQ GmbH
Konrad-Zuse-Ring 23
68163 Mannheim, Germany
www.1stq.de/en

