


611HPS Basis K hydrophilic



Technical specifications

Haigis-constants*	a0 = 1.839 a1 = 0.40 a2 = 0.10
A-constant*	118.9 (SRK/T)
Refractive Index	1.46
Abbe No.	58

*Stand IOL Con 12.2021

Lens type	One piece capsular bag IOL
Material	Biocompatible hydrophilic acrylic (25%)
Haptic	C haptic, 9° angulation
Optic	Spheric optic, 360° cell barrier
Filter	UV filter
Sterilisation	Steam
Injection system	1stInject 2.0HB, 1stInject 2.4HB

■ Basis function

Optical Power



611HPS Basis K hydrophilic

Article No.	Powers	Steps
BASW00	0.0 D to +9.0 D	1.0 D
BASW00	+10.0 D to +30.0 D	0.5 D



For the latest version,
scan the QR code.
www.1stq.de/basisk1074

**Do you have any
questions?**
export@1stq.eu



1stQ GmbH
Konrad-Zuse-Ring 23
68163 Mannheim, Germany
www.1stq.de/en