


Basis Q hydrophilic



Technical specifications

Haigis-constants*	$a_0 = -1.5646$ $a_1 = 0.3811$ $a_2 = 0.2265$
A-constant*	118.965 (SRK/T)
Refractive Index	1.46
Abbe No.	58

*IOL Con 12.2021

Lens type	One piece capsular bag IOL
Material	Biocompatible hydrophilic acrylic (25%)
Haptic	Q haptic, 0° angulation
Optic	Aspheric optic, 360° cell barrier
Filter	UV filter, blue light filter (optional)
Sterilisation	Steam
Injection system	1stInjct 1.8 to 30.0 D 1stInjct 2.0HB, 1stInjct 2.4HB

Basis function

Optical Power



Basis Q hydrophilic

Art. No. - clear	Art. No. - yellow	Powers	Steps
B2AW00	B2AWY0	0.0 D to +9.0 D	1.0 D
B2AW00	B2AWY0	+10.0 D to +30.0 D	0.5 D
S2AW00	S2AWY0	+31.0 D to +35.0 D	1.0 D



For the latest version,
scan the QR code.
1stq.de/basisq1072

**Do you have any
questions?**
export@1stq.eu



1stQ GmbH
Konrad-Zuse-Ring 23
68163 Mannheim, Germany
www.1stq.de/en